

Phonak Government Services Selection Guide | November 2026.

Hearing aids

RIC

Phonak Audéo™ E90




Phonak Audéo L90




Fitting ranges

S Receiver




Mild to moderately-severe hearing loss

M/M Sensor Receiver




Mild to severe hearing loss

P Receiver



Mild to severe hearing loss


UP Receiver



Moderate to profound hearing loss

Custom

Phonak Virto™ I90




*Wireless connectivity

Phonak Virto P90




All styles
40/109 (M)




For mild to moderate hearing loss

All styles
50/115 (P)




For mild to moderately-severe hearing loss

All styles
60/119 (SP)



For moderate to severe hearing loss


FS, HS, CA (312 battery)
70/127 (UP)




For moderate to profound hearing loss

BTE


Phonak Naída™ L90



Mild to moderately-severe hearing loss, all audiometric configurations




Mild to severe hearing loss, all audiometric configurations



Severe to profound hearing loss, all audiometric configurations

Phonak Naída Link M



Moderate to profound hearing loss, all audiometric configurations

Phonak Government Services Selection Guide | November 2026.

CROS



CROS E-R
Compatible with Audéo E90-R



CROS L-R
Compatible with Audéo L90-R/RT
and Naída L90-PR

Wireless Accessories



TV Connector 2.0



Phonak
PartnerMic™



Phonak
RemoteControl

Roger On™ V3

All-in-one microphone
and media streaming device



Colors

EON/Lumity



H0 Beige P1 Sand Beige P3* Sandalwood P4 Chestnut P5 Champagne P6 Silver Gray P7 Graphite Gray P8 Velvet Black 8A* Sienna Red

*N/A for Audéo or CROS Lumity

Naída Link M



P1 Sand Beige P4 Chestnut P6 Silver Gray P8 Velvet Black

Custom shells



26 Pink 14 Tan 22 Cocoa 28 Brown 06 Black R Red Transparent B Blue Transparent 21 Transparent

Custom faceplates



26 Pink 14** Tan 22 Cocoa 28** Brown 06 Black

** N/A on Virto Titanium

Virto Titanium faceplates†



26 Pink 22 Cocoa 06 Black

†Shell color U0 Titanium Gray

Earmolds⁺⁺/Domes



Phonak CROS Phonak Skeleton Phonak Semi-Skeleton Phonak Canal Phonak Helix Lock Phonak Canal Lock Phonak Half Shell Phonak Full Shell Carved Phonak Standard Full Shell

⁺⁺ Earmolds come in both silicone or acrylic



Cap Dome Open Dome Vented Dome Power Dome EasyGuard™ Open EasyGuard™ Vented EasyGuard™ Power SlimTip (Acrylic) SlimTip (Silicone) cShell cShell Titanium MAV TITNM SlimTip